

SIL-KRAS(G12C) protein standard – Arg/Lys ¹³ C ¹⁵ N labelling	
Human K-RAS (Uniprot accession number P01116-2 (1-169)) G12C mutant	
Product Description	Recombinant human K-RAS (1-169 sequence – G12C mutant) was expressed in <i>E. coli</i> with polyHis purification tag and purified using proprietary method including affinity chromatography.
Protein sequence	MSGSHHHHHHSSGIEGRMTEYKLVVVGACGVGKSALTIQLIQNHVDEYDPTI EDSYRKQVVIDGETCLLDILDITAGQEEYSAMRDQYMRTGEGFLCVFAINNTKS FEDIHHYREQIKRVKDSQVPMVLVGNKCDLPSRTVDTKQAQDLARSYGIPFIE TSAKTRQGVDDAFYTLVREIRKHKEK
Predicted Molecular Mass	21.16 kDa
Expression system	<i>E. coli</i>
Purity	Greater than 90% as determined by SDS-PAGE analysis
Labelling	ARG- ¹³ C ₆ , ¹⁵ N ₄ and LYS- ¹³ C ₆ , ¹⁵ N ₄ >99% as determined by LC-MS/MS analysis of trypsin digested SIL-protein
Tag information	PolyHis tag at the N-terminus
Protein Content	Quantitation is carried out by BCA protein assay
Formulation	Lyophilized from 20 mM HEPES, pH7.5, 150 mM NaCl and 1 mM DTT buffer. Reconstitute by adding ultrapure H ₂ O.
Shipping and storage	The product is supplied as a lyophilized powder and shipped at room temperature. Store at -80°C upon receipt. After reconstitution, protein can be kept at 4°C for a few weeks. Avoid repeated freeze/thaw cycles.
Usage	The product is for research use only. Not for diagnostic or therapeutic use.
Legal	This product is licensed under U.S. Patent No. 7,396,688 and foreign counterparts from E. I. du Pont de Nemours and Company. The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product for research and development only. Information about licenses for other uses is available from: E. I. du Pont de Nemours and Company; Attn: Associate Director, Commercial Development; DuPont Experimental Station E268; 200 Powdermill Rd.; Wilmington, DE 19803; 1-877-881-9787 (voice), 1-302-695-1437 (fax), licensing@dupont.com.

